

Jain Laboratory Instruments







Product Code. JL-TOMLE-10013

Cariollis Component of Acceleration

Description

Cariollis Component of Acceleration

This set up is designed to study coriollis component of Acceleration of a slider crank Mechanism.

The apparatus uses hydraulic analogy to represent the rotating slider.

It consists of a rotating block with two arms in opposite direction.

These tubes can be rotated at various speeds by using a swinging field motor, which also acts as a dynamometer to measure torque applied to rotating tubes.

A Perspex window on top cover helps to visualize the process. Rotameter is used to measure water flow rate through tubes.

Water is circulated by small monoblaock pump.

Specifications:-

Main Tank with fiberglass lining.

Rotating Arms 9/6 mm dia, 300 mm long.

Motor Swinging field, D.C., 0.5 HP.

Rotameter

Monoblock Pump

Control Panel comprising of:-
Speed Control Unit.
Speed Indicator.
Necessary switches.
Rigid support frame.
Range of Experiment:-
Coriollis Component of Acceleration can be determined at various speeds of rotation and water flow rates.
Service Required:-
A.C. Single Phase .230 V. stabilized supply.
Floor Space:-1.5 m X 1.5 m.
Contact JLab Export for your Educational School Science Lab Equipments. We are best educational equipment manufacturers india, educational equipments, educational equipments exporter, educational equipments exporters, educational equipments for science labs, educational lab equipment export.