



Product Code . JL-OI-10363

Perfect Tibial Nails

Description

Perfect Tibial Nails

Use 4.35mm Locking Screw for Nails Dia. 8mm & 9mm and

Use 4.8mm Locking Screw for Nails Dia. 10mm & 11mm and

Use 4.8mm Cancellous Locking Screw (Dual Core) for 3 most Proxima

Advanced Nail Design : New Anatomic bend for ease of nail insertion. Titanium Alloy for improved mechanical and fatigue properties.

Our Perfect Tibial Nail is indicated for fractures in the Tibial Shaft as well as for Metaphyseal and certain intraarticular fracture of the tibial head and pilon tibiae with :

Advanced Distal Locking Options : Distal Oblique Locking Option to prevent soft tissue damage and increase stability of the distal fragment.

Versatile Proximal Locking Option : Three unique and innovative locking options, in combination with Cancellous Dual Core locking Screws to increase the stability of the proximal fragment for proximal third fractures.

Medio Lateral Locking options enable primary compression or secondary controlled dynamization.

Same Nail for Right & Left Tibia

Use with Instrument Set - AS075I Holes.

Note: Distal Jig will work upto Nail Length 375 and beyond that use C-Arm both on the Nails Page and on the Jig Photos
