

Jain Laboratory Instruments







Product Code . JL-PALE-10529

Smoking burner and stand

Description

Smoking burner and stand

Description:-

The burner is put on and drum is smoked uniformly by rotating on the flame and a thin and black smoking is obtained. Then the drum is placed on the horizontal arm of smoking stand. The smoking stand enables the drum to be smoked on its own stand. A piece of glazed paper is pasted on the drum.

Smoking outfit comprises of basically 2 parts-Smoking stand & Smoking Burner kerosene or gas type.

Smoking Stand

Standardised to take any size of cylinders(6? or 10?) with stainless steel spindle, free moving in ball bearing with turning handle.

Smoking Burner

Kerosene type Smoking burner Gas Type.

