

Product Code . JL-PSAT-1209



Power Supply for Tetrode and Pentode Valve Model

Description

Power Supply for Tetrode and Pentode Valve Model

- Dual Power Supply for Pentode Valve Experiment, both stabiles $\pm 1\%$ workable on 220 Volt 50 Hz.
- With five meters.
- HT1 = 0 to 300 Volts 50 mA
- HT 2 = 0 to 300 Volts 25 mA
- Bias = 0 to 25 volts 5 mA
- L.T. = 6.3 Volt 3 Amps.
- Suitable for pentode/tetrode valve experiments.

Contact JLab Export for your Educational School Science Lab Equipments. We are best school science equipments suppliers, school science kits, school science lab equipment exporter, school science lab equipment manufacturer, school science lab equipment supplier, school science lab equipment suppliers.
