

## Jain Laboratory Instruments





Product Code . JL-S-2070



## **Reciprocating Shaking Machine**

## **Description**

Reciprocating Shaking Machine

Designed for shaking khan racks, bottles and flasks of different sizes by adjusting compartments.

Fitted with heavy duty motor to work on 220 V AC 50 Hz V A single phase.

The rate of motion is Adjustable mechanically by belt tension.

Speed is fixed in same models to approximately 275 or variable as detailed below.

Optional:- 0-60 minutes timer.

Size of Carrier Stokes/min

20"x13"x4" 250

31"x17"x4" 50-250

36"x18"x4" 50-250

Reciprocating Shaking Machine Manufacturer, Reciprocating Shaking Machine Suppliers, Reciprocating Shaking Machine India, Reciprocating Shaking Machine Exporter, buy Reciprocating Shaking Machine, Shaker, Reciprocating Shaking Machine, Analytical Research Equipments, buy Reciprocating Shaking Machine Online India

{ "@context": "http://schema.org/", "@type": "Product", "name": "Reciprocating Shaking Machine", "image": "https://www.jlabexport.com/images/catalog/product/1213005189\_2017-10-25.jpg", "description": "Reciprocating Shaking Machine Designed for shaking khan racks, bottles and flasks of

different sizes by adjusting compartments. Fitted with heavy duty motor to work on 220 V AC 50 Hz V A single phase. The rate of motion is Adjustable mechanically by belt tension. Speed is fixed in same models to approximately 275 or variable as detailed below. Optional:- 0- 60 minutes timer. Size of Carrier Stokes/min Reciprocating Shaking Machine Manufacturer, Reciprocating Shaking Machine Suppliers, Reciprocating Shaking Machine India, Reciprocating Shaking Machine Exporter, buy Reciprocating Shaking Machine, Shaker, Reciprocating Shaking Machine, Analytical Research Equipments, buy Reciprocating Shaking Machine Online India", "brand": "JLab Export", "sku": "5", "gtin8": "5", "gtin13": "5", "gtin14": "5", "mpn": "5", "aggregateRating": { "@type": "AggregateRating", "ratingValue": "5", "bestRating": "5", "worstRating": "0", "ratingCount": "5" } }