

Jain Laboratory Instruments





Product Code . JL-MI-2127

Digital Stop Clock

Description

Digital Stop Clock

To overcome the difficulties experienced in mechanical stop clocks.

JLab has introduced Quartz controlled 31/2 digit stop clock Start/Stop.

The operation is by means of a mini toggles witch Reset and by a push button Reset is possible only in Stop mode to avoid accidental resetting while the clock is running.

Range	999.90sec.
Resolution	0.1 sec.
Accuracy	±0.1%(Quartz controlled)
Display	12.5mm Bright seven segment display
Working Voltage	230 volts A.C. ±10%
Reset	By Push Button
Least Count	1 Second

Contact JLab Export for your Educational School Science Lab Equipments. We are best science educational equipments, science educational instruments, science lab equipments manufacturer, science laboratory instruments exporter, scientific educational equipments, scientific educational lab equipment.