

Product Code . JL-RE-2225

## Flocculator Jar Test Apparatus



### Description

---

Flocculator Jar Test Apparatus

#### Features:-

Accepts two, four & six samples

Digital indication of speed and time

Two pre-set programs

Built in illumination for easy viewing of samples

Flocculation/Jar Test Apparatus, is a flocculation (jar) tests on water and effluent samples.

A digital display clearly indicates the speed of the rotational stainless steel paddles, which can be varied from 25 to 250rpm.

The digital timer can be set to count down from 1 to 99 minutes.

After count down, the paddles stop and an alert sounds.

The timer then counts up in minutes in order to measure settling time.

Two pre-set programs allow routine speed and time parameters to be input and then selected by simply pressing one button.

Parameters are stored in the memory and can be recalled or altered at any time.

---

The intuitive touch sensitive keypad is totally flat for easy cleaning.

**Flocculator Jar Test Apparatus Manufacturer, Flocculator Jar Test Apparatus Suppliers, Flocculator Jar Test Apparatus India, Flocculator Jar Test Apparatus Exporter, buy Flocculator Jar Test Apparatus, Flocculator Jar Test Apparatus, Analytical Research Equipment, buy Flocculator Jar Test Apparatus Online India**

Specifications:-

No. of Jar	2	4	6
Speed range	25 to 200rpm		
Pre-set programs	2		
Timer	0 to 99minutes		
Controls Keypad	Touch sensitive		
Digital Display	LED		
Power	230V, 50Hz, 100W		

```
{ "@context": "http://schema.org/", "@type": "Product", "name": "Flocculator Jar Test Apparatus",  
"image": "https://www.jlabexport.com/images/catalog/product/369392829flocculator-jar-test-  
apparatus.jpg", "description": "Accepts two, four & six samples Digital indication of speed and time  
Two pre-set programs Built in illumination for easy viewing of samples Flocculation/Jar Test  
Apparatus, is a flocculation (jar) tests on water and effluent samples. A digital display clearly  
indicates the speed of the rotational stainless steel paddles, which can be varied from 25 to 250rpm.  
The digital timer can be set to count down from 1 to 99 minutes. After count down, the paddles stop  
and an alert sounds. The timer then counts up in minutes in order to measure settling time. Two pre-  
set programs allow routine speed and time parameters to be input and then selected by simply  
pressing one button. Parameters are stored in the memory and can be recalled or altered at any time.  
The intuitive touch sensitive keypad is totally flat for easy cleaning. Flocculator Jar Test Apparatus  
Manufacturer, Flocculator Jar Test Apparatus Suppliers, Flocculator Jar Test Apparatus India,  
Flocculator Jar Test Apparatus Exporter, buy Flocculator Jar Test Apparatus, Flocculator Jar Test  
Apparatus, Analytical Research Equipment, buy Flocculator Jar Test Apparatus Online India",  
"brand": "JLab Export", "sku": "5", "gtin8": "5", "gtin13": "5", "gtin14": "5", "mpn": "5",  
"aggregateRating": { "@type": "AggregateRating", "ratingValue": "5", "bestRating": "5", "worstRating":  
"0", "ratingCount": "5" } }
```

---