



Product Code . JL-NTM-224

Female Catheterization Training Simulator

Description

Female Catheterization Training Simulator

- The genital inserts are placed into the anatomically correct lower abdominal model and held in place with magnets.
 - The material is soft and flexible: the labia are movable and can be separated. This way the student can practice all the necessary gestures for catheterization (e.g. disinfection).
 - The material gives a realistic feeling when inserting and removing the catheter into the bladder with realistic resistance and insertion depths.
 - General urethral catheterization process can be operated.
 - There are resistance when advantage the catheter through urethral sphincter.
-