

Product Code . JL-BTD-3335

Electronic Micro pipettes



Description

Electronic Micropipettes

User friendly software.

Easy, user friendly calibration.

Recharging during operation.

Optimized ergonomics and working comfort.

Large liquid crystal display (LCD).

Built-in tip ejector.

Six models cover the complete pipetting range from 0.2 µl to 1000 µl.

Features

Model (Easy selection of modes from a single key)

Standard mode

Stepper mode

Dilution mode

Calibration mode

Variable working speeds

Speed can be selected on a scale of 1-5.

Volume Selection

Easy volume selection by 'UP' & 'DOWN' keys.

Large Fire Button

Single 'FIRE BUTTON' for filling and dispensing.

Large Liquid Crystal Display

Two Line Alphanumeric LCD.

Built-in streamlined tip ejector

For taken easy tip ejection and access to narrow naked bottles

Calibration

User friendly calibration mode facilitates easy calibration

Battery

9V NIHM rechargeable batteries are used.

8 hour of operation with fully charged **batteries**.
