

## ENT Doctor Chair

The function of up, down, lie, rise, and chair rotation(rear board angle with 90~180degrees) will be done by full automation mechanism

Arm rest can be rotated 360 degrees, giving the conformability to patients when getting on and down

Chair can be controlled by the switch in E.N.T. Unit

Head rest made from stainless steel mold are so solid and have semi-permanent life time without shaking , sinking & erosion

Seat cured by high frequency are so strong, endurable from scratch, not sewaged easily, which makes the patient's stability

Foot switch equipped basically has the function to operate the chair simply

Can be rotated 360 degrees by automation system

Long stroke supports the convenient treatment for children and adults with the best condition

Oil-pressure system adopted has semi-permanent life time