



**Product Code . JL-MS-3760**

## **Vertical Sterilizer**

### **Description**

#### **Vertical Sterilizer**

Corrosion & abrasion resistance

Excellent performance

Easy handling

#### **Specifications:**

Steam Working Pressure	1.26 kgf/cm <sup>2</sup> (2.2 kgf/cm <sup>2</sup> in case High speed sterilizer)
Steam Working Temperature	121oC (134o in case of High speed sterilizer)
Hyd. Tested at	Jacket-Twice the working pressure. Chamber-One & half time the working pressure.
Operating Voltage	400/440, 3 phase, AC supply, 50c/s.
Sterilization Period	20 to 25 minutes(in case of High speed sterilizer 5 to 7 minutes).
Steam Exhaust	5 to 7 minutes(in case of High speed sterilizer 1 minute)

---

Direct Contact Details  +91-8569909696  sales@jlabexport.com  
 [www.jlabexport.com](http://www.jlabexport.com)