

Product Code . JL-IP-3861

Syringe Infusion Pump



Description

Syringe Infusion Pump

Self-explanatory with abbreviated operating instructions on LCD.

Easy-to-read liquid-crystal displays all essential data in large size numbers.

Reducing the chance of erroneous readings.

LCD provided with back light for viewing in dark environments.

Selectable three-step occlusion detection pressures.

Low, mid & high-occlusion setting.

"Keep-Vein-Open" function.

Infusion rate in KVO is 1% but not less than 0.1 ml/hr.

This prevent the infusion line from blocking due to blood clots.
