

Product Code . JL-MV-3957



Respironics Lifecare Ventilator

Description

Respironics Lifecare Ventilator

Mode Selection: Control, Control + Sigh, Assist/Control, Assist/Control + Sigh, SIMV.

Tidal Volume: 0.05 to 0.20 \pm 0.02 liters; 0.20 to 3.00 liters \pm 10%.

Sign Volume: 150% of Set Tidal Volume.

Sign Frequency: 1 Sign Breath every 100 Tidal.

Breath Rate: 2 to 35 BPM \pm 0.5; 36 to 40 BPM \pm 2..

Aspiratory Flow Rate (Peak flow rate during aspiratory time): 10 to 120 LPM.

Battery Test Switch: Test Internal/External Battery.

Main Power Switch: On/Off, Recharge.

DC Circuit Breaker: Push/Pull Actuation.

AC Circuit Breaker: Push/Pull Actuation.

Airway Pressure Limit: 5 to 100 cm H₂O \pm 5.

Sensitivity: +18 cm H₂O \pm 1 to at least -6 cm H₂O.

Low Pressure Alarm: 2 to 40 cm H₂O \pm 2; 41 to 50 cm H₂O \pm 5.

Manual Sigh: Push to Activate One Breath.
