

Jain Laboratory Instruments







Product Code . JL-PLE-4260

Microtome Accessories

Description

Microtome Accessories

JLab pecial hone tone for microtome knife harpening (green) 10 pecial hone tone for microtome knife harpening (yellow) ize 10 leather trop made of best quality leather mounted on wooden frame, for final honing of microtome knives.

(Flock type) Embedding ring pencers tissue embedding rings have sturdy plastic mouldings which forms the permanent base of a wax block.

It is reinforced with dual bar center support to prevent warping these rings have a cross-hatch surface design that provide excellent paraffin bonding for clamping used base moulds including our plastic, metal aluminium disposable base moulds.

Embedding mould pencer permanent reusable metal base moulds for embedding cassettes are made of stainless steel and designed to ensure easy handling.

There are two different shapes available for various specimen sizes.

Embedding mould (disposable) pencers disposable moulds is designed to fit all existing cassettes.

Also as our base moulds are disposable you save time on cleaning and storing and they offer you less exposure to hazardous chemicals.

Et of three object holders for rotary microtomes.

Et of three different object holders for rocking microtomes.

| Razor attachment : a very simple attachment permitting a safety blade in place of microtome kni very suitable for cutting accurate sections. | |
|--|--|
| Micro abrasive suspension for coarse honing of microtome knives (500 gms.) | |
| Micro abrasive suspension for fine honing of microtome knives (500 gms.). | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |
| | |