



Product Code . JL-MI-4855

Capillary Tube Stand

Description

Capillary Tube Stand

Two types of tube stands are available i.e. 3 tubes and 6 tubes.

Then it is released to allow it to tighten against the corresponding tab on top of the stand.

In addition to this, tubes in our stands are positioned according to numbers for easy identifications.

Our stands are widely used to hold the micro sedimentation tubes by placing in the tube on the spring-loaded tab on the base.

We provide premium quality Capillary tube stands to our clients manufactured using sheet metal and finished in white enamel.
