

Product Code . JL-ELA-5395

Pipette Filler Pump



Description

Pipette Filler Pump

The thumb wheel can be rotated to precisely draw up or dispense liquids, or alternatively, the releaser may be pushed in for virtually instant emptying.

Acids and alkalis should be resisted.

It can be opened and cleaned easily.

Contact JLab Export for your Educational School Science Lab Equipments. We are best school science equipments suppliers, school science kits, school science lab equipment exporter, school science lab equipment manufacturer, school science lab equipment supplier, school science lab equipment suppliers.
