

## Jain Laboratory Instruments





Product Code . JL-ELA-5399

**Polishing Cloth** 

## Description

Polishing Cloth

A selection of premium polishing cloths for coarse and fine polishing of hard, soft and composite materials.

Most clothes are available with plain or PSA backing for easy mounting on the polishing discs.

The plain polishing cloths need to be clamped on the polishing discs with a clamping ring.

To facilitate this the plain discs are aproximately 2" larger than the polishing discs so the cloth folds over the edge of the disc for clamping.

8" plain polishing cloths are 10" diameter and 10" plain polishing cloths are 12" diameter.

The polishing cloths can be used with diamond compound, diamond spray or alumina suspension.

Hard Cloths

Low resilience, (e.g. Woven Nylon)

Used for rough and intermediate polishing

Soft Cloths

High resilience

Low nap (e.g. Rayon pile)

Used for final polishing

Contact JLab Export for your Educational School Science Lab Equipments. We are best school lab equipments supplier, school lab export, school lab exporter, school laboratory equipment manufacturers, school laboratory equipment suppliers, school science equipments.