

## Jain Laboratory Instruments







Product Code . JL-ELA-5474

## **Dry Cleaning Apparatus**

## **Description**

**Dry Cleaning Apparatus** 

Dry cleaning Apparatus is used for preparing fabric samples for flammability testing.

This apparatus can also be used for testing color fastness and dimensional changes during dry cleaning and during laundering.

Dry Cleaning Apparatus consists of a motor driven corrosion resistant stainless steel tumbler, which rotates around an axis 50anddeg; inclined to tumbler axis at approximately 45 rpm.

The tumbler is also provided with a stainless steel cover with a rubber seal and a quick acting positive clamp.

The tumbler is rotated with the help of an electric motor through pulley-belt, reduction gearbox, and chain and sprocket drive.

Contact JLab Export for your Educational School Science Lab Equipments. We are best school lab equipments supplier, school lab export, school lab exporter, school laboratory equipment manufacturers, school laboratory equipment suppliers, school science equipments.