

Product Code . JL-E-550

Brownlee Electrolysis Apparatus



Description

Brownlee Electrolysis Apparatus

The most popular form of Electrolysis Apparatus.

Consists of Electrodes attached to insulated connecting rods which are in-turn connected to terminals.

Mounted on a non-conducting plastic holder, designed to rest on top of a glass jar / beaker.

The plastic holder has clips to hold two test tubes.

Operates on 6 - 10v DC, power supply not included.

Supplied complete with:-

Beaker (1000mL)

Graduated Test Tubes

Un-graduated Test Tubes

Rubber Stoppers for Test Tubes

A preview of the included instructions can be found at the bottom of this listing.

Measurements are approximate.

Instructions Included.

Shipping Weight 5.00 lbs
