

Description

AIRWAY SIMULATING:

- Compression indicated by indicating lights.
- Counter and audio instruction.
- Which is located in the central position in the chest, with green flash.
- Compression times is indicated by the counter.
- Which is located on the left side of the monitor.
- When the compression position is right on the CPR Manikin.

SOME PROCEDURES WILL BE GUIDED BY AUDIO INSTRUCTION:

- The ratio between compression and mouth-to-mouth ventilation on this CPR Manikin is: 30:2
 - the numbers of a whole circle: five circles of 30 times compression and 2 times ventilation.
 - Report: A report is printed after practice.
 - Mode: Training mode; Test mode.
 - Pupil observation: From the beginning of the practice to the end, people can observe the change of the pupil vary from release to collect.
 - To check the pulse of the cervical artery. # Power supply: 220V AC puts in through power adapter, and the output is 24V DC. The power supply for the Report Printer is 9V.
 - Email us at sales@jaincoexport.com to further inquire about the product.
 - Visit www.jaincoexport.com for further details.
 - Also you can reach us at +91-85699-09696 (Whatsapp).
-