

Jain Laboratory Instruments







Product Code. JL-APE-7308

Tel tron Tube Holder S

Description

Tel tron Tube Holder S

The Auxiliary Coil can be inserted into the stepped front side of the tube holder.

An easily replaceable fuse is located in the holder to protect the valuable tubes from excessive voltage.

The slots for holding coils have metric rulings to ensure precise adjustments with a maximum spacing of 150 mm.

This holder is the mounting point for the entire line of Tel tron Student ('S') Series and related accessories, such as the Helmholtz Coils S and Auxiliary Coil.