

Jain Laboratory Instruments







Product Code . JL-H-746

Ether Thermoscope Pulse

Description

Ether Thermoscope Pulse

- Glass construction.
- On holding the lower bulb in hand causes the liquid ether inside to boil and shoot into the top bulb.
- The body heat causes the liquid to appear as if boiling.
- This is a very fragile and delicate product.
- Any excess amount of pressure on the product can cause it to shatter or break.
- Use extreme caution.
- DO NOT IN ANY CASE HEAT THIS PRODUCT, IT WILL CAUSE IT TO EXPLODE.
- The liquid inside is extremely volatile and inflammable.
- We are not liable for any kind of injury, loss or damage of any kind, to property or to you or anything or anyone.
- When you buy this product you expressly agree to this statement.
- A preview of the included instructions can be found at the bottom of this listing.
- Measurements are approximate.