

## Jain Laboratory Instruments





## Product Code . JL-SI-8143

## **Groove Type Rheostat Tube**

## Description

Groove Type Rheostat Tube

The special brass slider moves the full length of the Rheostat, moves rapid with ease.

Contact with Eureka resistance wire is made by two copper brushes having large contact area to give smooth current control and very little wear.

The only difference is that imported bare with Eureka wire (having high stability and strength) is wound on a thick wall porcelain/cement pipe having a deeply grooved outer surface to provide air conducts for convection cooling.

Contact JLab Export for your Educational School Science Lab Equipments. We are best physics lab equipment supplier, physics lab equipment suppliers in ambala, physics lab equipments manufacturers, physics lab equipments suppliers, physics labware, school educational equipments.