

Product Code . JL-SM-8863

Hand Operated Casagrande Equipment



Description

Hand Operated Casagrande Equipment

Use:-

Measuring the Liquid limit.

Specification:-

Hand operated on the counter.

A brass cup suspended from a carriage designed to control its drop exactly 1cm onto a hard rubber base.

Adjustable height of drop by a nut at the rear of the cup.

Casagrande grooving tool.

Rubber feet isolate base from work surface.

Contact JLab Export for your Educational School Science Lab Equipments. We are best scientific laboratory equipments manufacturer, scientific laboratory equipments supplier, technical educational equipment exporter, technical educational equipment manufacturer, technical educational equipment supplier, technical educational equipments exporter.
