

**Product Code . JL-SM-8866**

## Cone Apparatus



### Description

---

Cone Apparatus

Use:-

For accurately determining the density of gravels and coarse soils.

Specification:-

305mm (12") diameter cone apparatus.

The unit has two identical cones with the valve in the centre.

Metal cylinder attached to the top cone with a clear plastic end for viewing sand flow and handles for carrying.

Accessories:-

Circular density plate.

Poly sample bags, qty. 100.

---

25°-125°F pocket dial thermometer.

Bristle brush.

1 in steel chisel.

Stainless steel spoon.

Rubber Mallet.

Contact JLab Export for your Educational School Science Lab Equipments. We are best technical educational equipments manufacturers, technical educational equipments manufacturer, technical educational equipments supplier, technical educational instruments, technical educational instruments exporter, technical educational instruments manufacturer, vocational training lab instruments exporter.

---