

Jain Laboratory Instruments







Product Code . JL-SM-8887

Triaxial Test Set

Description

Triaxial Test Set

Use:-

For the Tri-axial test; Unconsolidated Undrained (UU).

Consolidated Undrained (CU).

Consolidated Drained (CD).

Unconfined Compression (UC).

Triaxial Cell:-

Triaxial cell for specimens up to 70 mm diameter.

Pedestal for 70 mm specimens.

Top cap with drainage for specimens of 70 mm.

The base plate, in aluminum for specimens of 70 mm.

Porous disc for 70 mm specimens (set of 2).

Suction sleeve stretcher for 70 mm - 2.8" - 3" specimens. Membrane, rubber latex for 70 mm - 2.8" specimens. O-ring neoprene for 70 mm - 2.8" specimens. Filter drain for 70 mm specimens pack of 50. Two-part split mould for 70 mm specimens. Contact JLab Export for your Educational School Science Lab Equipments. We are best school science lab equipment suppliers in india, school science lab equipments exporter, school science lab equipments manufacturers, school science lab equipments manufacturer, school science lab equipments supplier, school science lab equipments suppliers.