



Product Code . JL-AMLE-8941

Flywheel Apparatus

Description

Flywheel Apparatus

Self contained.

Set of weights.

Stop watch supplied.

Running in ball bearings.

A load hanger and cord supplied.

Instruction manual for student and lecturer provided.

To demonstrate the Second law of motion and energy storage.

Wall mounting includes a pointer to the periphery of the wheel.

250 mm x 30 mm steel flywheel with engraved line on periphery

Dimensions and Weight (excluding set of weights):-

215 x 280 x 250 mm ; 19 kg (approx)
