

Jain Laboratory Instruments







Product Code. JL-DOMLE-9141

Corollis Component of Acceleration

Description

Corollis Component Of Acceleration

Coriolli's components of Acceleration can be determined at various speeds of rotation at water flow rates.

Parts/Component Details:-

Rotameter 300 to 3000 LPH.

Main Tank Fabricated out of fiberglass plastic sheet.

Electrical Motor D C. Swinging field, 0.5 H.P. 50 RPM.

Rotating Arms 9 mm/6 mm orifice diameter, 300 mm long.

Mono-block Pump Single phase, pump with Motor 2400 LPM discharge.

Control Unit consisting of:-

Speed Control,

Speed Indicator ,			
Main Switch.			
Rigid Support Structures.			