

Jain Laboratory Instruments







Product Code . JL-MAI-9835

Dark Skin Cancer Trainer

Description

Dark Skin Cancer Trainer

Continued exposure to sun-light can cause damage to skin cells.

If the damaged cells do not die or repair themselves, they degenerate and visible skin cancer develops.

The Life/form Skin Cancer Trainer has been developed to assist with educating healthcare providers about the various types of skin cancer and what they may look like.

The enlarged trainer shows a normal mole, basal cell carcinoma (basalioma), squamous cell carcinoma (spinalioma) and malignant melanoma.

Variations of how the different types of skin cancer may appear are shown.

The Life/form Skin Cancer Trainer includes trainer, sturdy base, two-sided informational key card, easel and storage box.

Size:- 8-1/2 in. x 11 in. x 1/4 in. (1/2 in. with board).