

Jain Laboratory Instruments



Product Code . JL-PLE-9888



Worcester Current Balance Kit

Description

Worcester Current Balance Kit

An elastic band, zeroing is effected using wire riders.

The inducation is made on a vertical wood strip positioned by a wood block.

Each correct balance consists, when assembled, of a coil which, when conducting, acts on a small magnet attached to the end of a drinking straw pointer with a needle pivot on an aluminium channel support.

4 mm reel of 16 mm self adhesive cellulose tape, 1 ppick 20 needles, 4 wood blocks, 10 wood strips.

Kit comprises 50 drinking straws, supply of elastic bands, 8 alcomax magnets, 55 gm reel of tinned copper wire 0.45 mm diameter, 4 wood bases carrying coil, channel and terminals.

Complete with instructions.