



**Product Code . JL-SW-1305**

## **Resonance Apparatus Complete Advanced**

### **Description**

---

#### **Resonance Apparatus Complete Advanced**

- To investigate Resonance of Air Columns, calculate Velocity of Air, Wavelength of Tuning fork.
- Included are heavy Support Stand with rod, two Clamps.
- Graduated Acrylic tube, Aluminum reservoir, plastic tubing, Tuning Fork Set of 8 and one Thermometer Red Spirit Filled.
- Complete with detailed Instructions.

#### **Cast Iron Support Stand:**

- Base- 8.7 inch (220 mm) x 5.8 inch (145 mm), Rod- 36 inch (910 mm) Long x 0.5 inch (12.5 mm) dia.

#### **Graduated Acrylic Tube:**

- 39 inch (1000 mm) long x 1.5 inch (38 mm) dia.

#### **Graduations:**

- 0 to 900 mm x 1 mm

#### **Tuning Forks Set of 8:**

- 1 each of C-256, D-288, E-320, F-341.3, G-384, A-426.6, B-480, C-512.
- Length Varying between 4.5 inch (120 mm) to 6 inch (150 mm)

---

**Thermometer:**

- -10 to 110C, 12 inch (300 mm)

**Instruction Manual:**

- A preview of the included instructions can be found at the bottom of this listing.
-