

Product Code . JL-HBSR-1986

Semi Fowler bed (Super)



Description

Semi Fowler bed (Super)

- Two section top of uniformly perforated CRCA sheet.
- Backrest section maneuvered by the help of Actuators to achieve smooth and efficient operation without any jerks.
- ABS bows easily removable for access to the patient
- The bed is fitted with 125 mm die swiveling castors with center locking
- Frame Work made of rectangular M.S. tubes
- Perforated CRC sheet top in two sections
- Backrest positions obtained by electric actuator operator by hand remote system
- Telescopic IV Rod with Four Locations
- S.S collapsible safety side railings
- Bed mounted on 100mm deluxe castors, two with brakes
- ABS molded head and foot bows
- Pre Treated & Epoxy Powder Coated Finish

approximate Size

{2160(L)× 1015(W) ×600(H)} mm.

mattress size

{1980(L)×910(W)} mm.
