

Product Code . JL-FTA-2216

Digital Friability Test Apparatus



Description

Digital Friability Test Apparatus

Features:-

Seven segment LED display

Sturdy motor to drive drums at constant speed

Digital Friability Test Apparatus is a precision instrument suitable for measuring obression strength of tablets.

The unit is equipped with two transparent acrylic drums which rotate at a speed of 25 ± 1 rpm.

Each of the two acrylic drums are provided with an arm which carries the tablets along with it up to a predetermined height and allows them to fall from that specified height, while the drums are rotating,

It is fitted with a digital revolution counter showing the current number of revolutions for the test.

A thumb wheel switch is used to set the desired number of rotations.

Provision is made to remove the lid of the acrylic drum by unscrewing just one knob.

This is available in 1 drum and 2 drum assemblies.

Digital Friability Test Apparatus Manufacturer, Digital Friability Test Apparatus Suppliers, Digital Friability Test Apparatus India, Digital Friability Test Apparatus Exporter, Friability Test Apparatus, buy Digital Friability Test Apparatus, Digital Friability Test Apparatus, Analytical

Specifications:-

Display	3 digit bright red seven segment LED display
Accuracy	± 1 revolution
Resolution	1 revolution
Range	0-999 revolutions
Operating Temp. Range	10°-50° C
Power	230 V ± 10%, 50 Hz. AC
Dimension	230 mm (H) X 250 mm (W) X 540 mm (L) (Approx.)
Weight	9.5 Kg. (Approx
Accessories	Drum 1 No., 2 Nos., Dust Cover-1 No.

```
{ "@context": "http://schema.org/", "@type": "Product", "name": "Digital Friability Test Apparatus",  
"image": "https://www.jlabexport.com/images/catalog/product/1100669926digital-friability-test-  
apparatus.jpg", "description": "Seven segment LED display Sturdy motor to drive drums at constant  
speed Digital Friability Test Apparatus is a precision instrument suitable for measuring obression  
strength of tablets. The unit is equipped with two transparent acrylic drums which rotate at a speed of  
25 ± 1 rpm. Each of the two acrylic drums are provided with an arm which carries the tablets along  
with it up to a predetermined height and allows them to fall from that specified height, while the drums  
are rotating, It is fitted with a digital revolution counter showing the current number of revolutions for  
the test. A thumb wheel switch is used to set the desired number of rotations. Provision is made to  
remove the lid of the acrylic drum by unscrewing just one knob. This is available in 1 drum and 2  
drum assemblies. Digital Friability Test Apparatus Manufacturer, Digital Friability Test Apparatus  
Suppliers, Digital Friability Test Apparatus India, Digital Friability Test Apparatus Exporter, Friability  
Test Apparatus, buy Digital Friability Test Apparatus, Digital Friability Test Apparatus, Analytical  
Research Equipment, buy Digital Friability Test Apparatus Online India", "brand": "JLab Export",  
"sku": "5", "gtin8": "5", "gtin13": "5", "gtin14": "5", "mpn": "5", "aggregateRating": { "@type":  
"AggregateRating", "ratingValue": "5", "bestRating": "5", "worstRating": "0", "ratingCount": "5" } }
```
