

## Jain Laboratory Instruments







Product Code . JL-TMM-2451

## Tablet Making Machine Hand Operated

## **Description**

Tablet Making Machine (Hand operated)

Useful for making tablets up-to 100 mg.

Size of punch & die of 9 mm.

Tablet Making Machine Hand Operated Manufacturer, Tablet Making Machine Hand Operated Suppliers, Tablet Making Machine Hand Operated India, Tablet Making Machine Hand Operated Exporter, Tablet Making Machine Hand Operated, Analytical Research Equipments, buy Tablet Making Machine Hand Operated Online India

Automatic compression and ejection device has been provided, actuated by revolution of the wheel drive.

Feed and removal of tablets has also been made automatic.

Eccentric pressure system has been applied for easy setting of thickness and hardness of tablets.

A stainless steel hopper has been given.

The unit has been mounted on a mild steel fabrication.

Completed with one die and punch set made of chromium steel.

{ "@context": "http://schema.org/", "@type": "Product", "name": "Tablet Making Machine Hand Operated", "image": "https://www.jlabexport.com/images/catalog/product/359188680tablet-making-machine-hand-operated.jpg", "description": "Useful for making tablets up-to 100 mg. Size of punch & die of 9 mm. Tablet Making Machine Hand Operated Manufacturer, Tablet Making Machine Hand Operated Suppliers, Tablet Making Machine Hand Operated India, Tablet Making Machine Hand Operated Exporter, Tablet Making Machine Hand Operated, Analytical Research Equipments, buy Tablet Making Machine Hand Operated Online India Automatic compression and ejection device has been provided, actuated by revolution of the wheel drive. Feed and removal of tablets has also been made automatic. Eccentric pressure system has been applied for easy setting of thickness and hardness of tablets.", "brand": "JLab Export", "sku": "5", "gtin8": "5", "gtin13": "5", "gtin14": "5", "mpn": "5", "aggregateRating": { "@type": "AggregateRating", "ratingValue": "5", "bestRating": "5", "worstRating": "0", "ratingCount": "5" } }