



Product Code . JL-MS-3760

Vertical Sterilizer

Description

Vertical Sterilizer

Corrosion & abrasion resistance

Excellent performance

Easy handling

Specifications:

Steam Working Pressure	1.26 kgf/cm ² (2.2 kgf/cm ² in case High speed sterilizer)
Steam Working Temperature	121°C (134°C in case of High speed sterilizer)
Hyd. Tested at	Jacket-Twice the working pressure. Chamber-One & half time the working pressure.
Operating Voltage	400/440, 3 phase, AC supply, 50c/s.
Sterilization Period	20 to 25 minutes(in case of High speed sterilizer 5 to 7 minutes).
Steam Exhaust	5 to 7 minutes(in case of High speed sterilizer 1 minute)