



Product Code . JL-NP-4023

Radiant Warmer with Under Surface Phototherapy

Description

Radiant Warmer with Under Surface Phototherapy

Easy to use

Focused heating

Servo control mode

Air control mode

Manual mode

Under Surface Unit:

Detachable source unit can be used along with radiant warmer & infant trolley

Two, 2 feet TLD daylight tubes

Two, 2 feet imported low pressure mercury vapour fluorescent tubes emitting blue radiation with a peak at 450 nm specially used for hyperbilirubinaemia treatment (TL 20W/52)

M.S. fabrication with powder coated paint

Clear acrylic filter

Source cooling fan

