

Jain Laboratory Instruments





Product Code . JL-LOL-4024

LED OT Lights

Description

LED OT Lights

Radio graphic output: 40kw as per (Parts)

Tube filament current stabilization through constant voltage transformer (CVT)

Full wave rectification

Digital display of KV, mA and MAS

Overload protection (Blockage of Exposure)

Shortest exposure time 10 MAS in 500mA

Supports serialography up to 6fps kv, mA and MAS values set with step switches fluoroscopy time totaliser

5 LEO indicators for Line, 0/L exposure, ready and bucky in voltage stabilization (compensation)

Through course and fine voltage compensation switches voltmeter and mA meter for viewing input line voltage and ma during exposure