

Jain Laboratory Instruments







Product Code . JL-SL-4770

Clinostat Clock Type

Description

Clinostat Clock Type

Technically sound Clinostat Clock Type can be purchased here for laboratory demonstrations and experiments. Clinostat Clock Type Manufacturer India, Clinostat Clock Type Exporters, Physics Instruments Clinostat Clock Type, Clinostat Clock Type Suppliers, Buy Clinostat Clock Type, Buy Clinostat Clock Type India, Didactic, Didactic Physics Instruments, Physics Instruments Exporters India, Physics Instruments Suppliers, Physics Instruments Manufacturer. It is designed to take a revolution after every hour for 24 hours in both horizontal and vertical positions. The plant is seeded in the pot, with clockwise movement on densely crafted metal base.

Winding mechanism is driven by hands.

Contact JLab Export for your Educational School Science Lab Equipments. We are best educational laboratory instrument suppliers, educational laboratory instruments exporters, educational laboratory microscopes, educational laboratory suppliers, educational physics lab equipments, educational science equipments.

{ "@context": "http://schema.org/", "@type": "Product", "name": "Clinostat Clock Type", "image": "https://www.jlabexport.com/images/catalog/product/1993930291clinostat-clock-type.jpg", "description": "Clinostat Clock Type Technically sound Clinostat Clock Type can be purchased here for laboratory demonstrations and experiments. Clinostat Clock Type Manufacturer India, Clinostat Clock Type Exporters, Physics Instruments Clinostat Clock Type, Clinostat Clock Type Suppliers, Buy Clinostat Clock Type, Buy Clinostat Clock Type India, Didactic, Didactic Physics Instruments,

Physics Instruments Exporters India, Physics Instruments Suppliers, Physics Instruments Manufacturer. It is designed to take a revolution after every hour for 24 hours in both horizontal and vertical positions. The plant is seeded in the pot, with clockwise movement on densely crafted metal base. Winding mechanism is driven by hands.", "brand": "Jlab Export", "sku": "5", "gtin8": "5", "gtin13": "5", "gtin14": "5", "mpn": "5", "aggregateRating": { "@type": "AggregateRating", "ratingValue": "5", "bestRating": "5", "worstRating": "0", "ratingCount": "5" } }