

## Jain Laboratory Instruments





Product Code . JL-AE-6715

**Tulgren Funnel** 

## Description

Tulgren Funnel

A soil sample is placed on the perforated plate.

A 220V heating lamp is held above the soil, driving the organisms down through the funnel into a collecting jar.

Supplied with coarse and fine mesh and 3 core cable.

Requires a low wattage bulb and a beaker (not provided).

Retort Rod and Base not included.

Tulgren Funnel India, Tulgren Funnel Manufacturer, Biology Lab Equipments, Biology Equipment Tulgren Funnel, Tulgren Funnel Suppliers, Tulgren Funnel Export, Educational Biology Equipments, Animal Equipments Tulgren Funnel, Tulgren Funnel Online, buy Tulgren Funnel, Tulgren Funnel Manufacturer, buy Tulgren Funnel Online India.

Contact JLab Export for your Educational School Science Lab Equipments. We are best educational laboratory instrument suppliers, educational laboratory instruments exporters, educational laboratory microscopes, educational laboratory suppliers, educational physics lab equipments, educational science equipments.

{ "@context": "http://schema.org/", "@type": "Product", "name": "Tulgren Funnel", "image":

"https://www.jlabexport.com/images/catalog/product/116259677TulgrenFunnel.jpg", "description":

"Tulgren Funnel A soil sample is placed on the perforated plate. A 220V heating lamp is held above

the soil, driving the organisms down through the funnel into a collecting jar. Supplied with coarse and fine mesh and 3 core cable. Requires a low wattage bulb and a beaker (not provided). Retort Rod and Base not included. Tulgren Funnel India, Tulgren Funnel Manufacturer, Biology Lab Equipments, Biology Equipment Tulgren Funnel, Tulgren Funnel Suppliers, Tulgren Funnel Export, Educational Biology Equipments, Animal Equipments Tulgren Funnel, Tulgren Funnel Online, buy Tulgren Funnel, Tulgren Funnel, Tulgren Funnel Online India.", "brand": "Jlab Export", "sku": "5", "gtin8": "5", "gtin13": "5", "gtin14": "5", "mpn": "5", "aggregateRating": { "@type": "AggregateRating"; "0", "ratingCount": "5" } }