

**Product Code . JL-SI-8063**

## Callender & Griffith Bridge



### Description

---

Callender & Griffith Bridge

Jockey slides on brass rod.

It is an improvement over the potentiometer.

The 50 cm long slide wire of 1 ohm resistance is calibrated into 100 equal parts.

Three in one - Carry's Foster Bridge, Potentiometer & Callender & Griffith's Bridge.

14 coils are arranged in series in dials from each being of 1 ohm resistance terminals arranged for appropriate working of different instruments.

Contact JLab Export for your Educational School Science Lab Equipments. We are best science educational equipments, science educational instruments, science lab equipments manufacturer, science laboratory instruments exporter, scientific educational equipments, scientific educational lab equipment.

---