

Product Code . JL-CP-8460

Elastrator



Description




Elastator

9-1/2" L. 2" opening.

Opens up to 2" in diameter.

The specially compounded rubber ring is easily put in place with the Elastrator, and left there.

The circulation is shut off, and any sensation of discomfort passes quickly as the part becomes numb, and eventually falls off.

Jain Laboratory Instruments Pvt. Ltd,
Hargolal Road, Ambala Cantt, Haryana India
Direct Contact Details  +91-8569909696  sales@jlabexport.com
 www.jlabexport.com